

Chapter 22: Introduction to Electroanalytical Chemistry


Two general categories:

- 1) Potentiometric Systems – measure voltage (i.e., potential) of a galvanic cell (produces electricity spontaneously)
- 2) Voltammetric Systems – control potential & usually measure current in an electrolytic cell (consumes power to cause an electrochemical reaction to occur)


Potentiometry

- Determine concentrations by measuring the potential (i.e., voltage) of an electrochemical cell (**galvanic cell**)
- Two electrodes are required
 - 1) Indicator Electrode – potential responds to activity of species of interest
 - 2) Reference Electrode – chosen so that its potential is independent of solution composition.

Representation of Electrochemical Cell



Representation of Electrochemical Cell



Ⓟ - Represents device to measure potential (voltage) without drawing significant current i.e potentiometer or electrometer (high input impedance $\geq 100 \text{ M}\Omega$ (mega ohms))

$$E_{\text{cell}} = E_{\text{ind}} - E_{\text{ref}} (+ E_{\text{J}})$$

E_{J} = junction potential, a non-ideal potential which develops across the interface between two dissimilar solutions

$$E_{\text{cell}} = E_{\text{ind}} - E_{\text{ref}} (+ E_{\text{J}})$$

Nernst Equation

$$E = E^{\circ} - \frac{RT}{nF} \ln \frac{[\text{Red}]}{[\text{Ox}]}$$


Where

- R = gas constant
- T = absolute temperature
- n = number of electrons in reaction
- F = Faraday's constant
- E = potential
- E° = standard potential
- [Red] = molar concentration of reduced form of species
- [Ox] = molar concentration of oxidized form of species

Reference Electrodes

- The Normal Hydrogen Electrode (NHE) is important historically and could serve as a reference electrode today, however, it is impractical, requiring a source of H₂ gas at constant pressure, and is highly flammable.

The Hydrogen Electrode





- Hydrogen gas is bubbled over an inert platinum electrode



Reference Electrodes

- The Calomel Electrode or Saturated Calomel Electrode (SCE) is the next most important reference electrode historically and was used almost exclusively for many decades as the reference electrode of choice
- Calomel is the insoluble compound Hg_2Cl_2
- The electrode half reaction is
- $\text{Hg}_2\text{Cl}_2 + 2 e^- \leftrightarrow 2 \text{Hg} + 2\text{Cl}^- \quad E^\circ = 0.242 \text{ v}$

Calomel Reference Electrode



Calomel Reference Electrode

- Can use 1 M or 0.1 M KCl rather than a saturated solution
- E for reference changes slightly with any change in concentration from the Nernst Eq
- Temperature coefficient of reference electrode is less with 1 M or 0.1 M than for SCE
- SCE often gets clogged if solution dries out